TECHNOTE 106

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Procedure for EDC/NHS Crosslinking of Carboxylates with Primary Amines**A: Materials Required**

- ☐ Activation Buffer: 0.1 M MES (2-[morpholino]ethanesulfonic acid), 0.5 M NaCl, pH 6.0.
- ☐ Phosphate-buffered Saline (PBS): 0.1 M sodium phosphate, 0.15 M NaCl, pH 7.2-7.5
- ☐ Protein: Prepare 1 ml of Protein in activation buffer at ~10 mg/ml

- Hapten: dissolved at 1-10 mg/ml in PBS or other amine-free buffer, pH 7-8
- EDC (1-ethyl-3-[3-dimethylaminopropyl]carbodiimide)

B. NHS-ester Activation

1. Add 0.4 mg of EDC (final concentration 2 mM) directly to 1 ml of Protein, which based on a 50 kDa protein, results in a 10-fold molar excess of EDC to Protein.
2. Add either 0.6 mg of NHS or 1.1 mg of Sulfo-NHS to the reaction (final concentration 5 mM).
3. Mix reaction components well and react for 15 minutes at room temperature.
4. (Optional): Separate activated Protein from excess EDC, EDC-byproducts, NHS and (if used) 2-mercaptoethanol using an appropriate size desalting column that has been equilibrated with PBS. Follow desalting column instructions and recover the fraction containing the activated protein. If using absorbance at 280 nm to identify fractions containing protein, be aware that NHS and Sulfo-NHS absorb strongly at 260-280 nm.

C. Amine Reaction

1. If step B.4 was not performed (i.e., buffer not exchanged using a desalting column), then increase buffer pH above 7.0 using concentrated PBS or other non-amine buffer such as sodium bicarbonate.
2. Add Hapten to the solution containing activated Protein.
3. Mix the solution well and then allow reaction to proceed for 2 hours at room temperature.
4. (Optional): Quench reaction by adding hydroxylamine to a final concentration of 10 mM. The excess hydroxylamine reacts to all NHS esters remaining on the surface of Protein, resulting in conversion of the original carboxyl groups to a hydroxamic acid. Alternative quenching reagents include 20-50 mM Tris, lysine, glycine and ethanolamine. Addition of base to raise the pH > 8 will promote hydrolysis of the NHS esters, thereby regenerating the original carboxyl groups.

Schematic diagram